WILLIAMSOUND

PPA R37-8 72 MHz FM Receiver

USER MANUAL

Safety Information

HEARING SAFETY

CAUTION!

This product is designed to amplify sounds to a high volume. To protect your hearing and the hearing of others:

1. Make sure the unit is OFF or the volume is turned down before putting on the earphones.
2. Set volume level at the minimum setting you need to hear.
3. Do not allow children or other unauthorized persons to have access to this product.

BATTERY SAFETY

CAUTION!

Do not attempt to recharge alkaline or other non-rechargeable batteries, which may explode, release dangerous chemicals, cause burns, or cause other serious harm to the user or product.

MEDICAL DEVICE SAFETY

CAUTION!

Before using this product with a pacemaker or other implantable medical device, consult your physician or the manufacturer of the device to make sure that you are using this product in accordance with their safety guidelines.

Operating Instructions

4. Open battery compartment door on back of receiver by pressing the tab down.
5. If using rechargeable NiMH batteries, slide the Alkaline/ NiMH switch to the NiMH position (Figure 1). Default setting is Alkaline for non-rechargeable batteries.



Figure 1: Battery selection switch

6. Insert 2 AA batteries over the fabric strip. Note correct polarity (+/-).
7. Insert earphone or neckloop into the headphone jack (Figure 2).
8. Choose the channel by rotating the channel selector knob to the desired channel number (Figure 2).
9. Turn Receiver ON and adjust volume control as needed (Figure 2).

Figure 2: PPA R37-8 Controls

- The **“ON”/LOW BATTERY INDICATOR** (Figure 2) is used to indicate: power, low batteries or sleep mode. Refer to flash codes chart on the next page.
- The **stereo headphone jack** accepts either stereo or mono headphones or earphones.
- The **adjustable squelch** (Figure 1) can help reduce user fatigue by turning off the audio when the signal becomes weak. If the operator is annoyed or fatigued by noise, turn the squelch up. This will reduce the effective range of the receiver as a trade-off to listener comfort. Similarly if the operator needs more range and can tolerate the noise, turn the squelch down*.

To adjust the squelch: using a small screwdriver, turn the knob to the right for more squelch; left for less squelch.

*As the radio signal decreases, the noise present at the headphones increases. The squelch circuit senses the signal level and turns off the audio as the noise increases. More squelch means this trigger point occurs earlier; less squelch means it occurs later (or not at all, with no squelch). Apparent range is also affected. More squelch reduces the apparent range; less squelch increases the apparent range.

PPA R37-8 LED Flash Codes

LIGHT CODE	INDICATES	REMEDY
On-Off-On-Off (continuous)	Low Batteries	Replace or recharge batteries
Fade-on, Fade-off (continuous)*	No transmitter on current channel*	After 6 min, Sleep/ power save mode starts (extends battery life)*
On Steady	Unit is operating normally	n/a

*The R37-8 returns to normal operation automatically when a transmitter is re-established on the same channel.

Belt Clip Installation

For hands-free operation of the PPA R37-8 receiver, use the provided belt clip to secure the unit on your belt, pocket or pants.

To Install: Position the belt clip on the rear of the PPA R37-8 receiver as shown in Figure 3A. Turn the belt clip 180° left or right as shown in Figure 3B. The belt clip is now installed and ready for use.

To Remove: Turn the belt clip 180° so the edge points toward the top of the unit as shown in Figure 3A. Gently pull the belt clip away from the unit to remove.

Figure 3: Belt Clip Installation


Figure 3A


Figure 3B

System Specifications

PPA R37-8 Receiver

Dimensions:	4.5" H x 2.85" W x 1.38" D (115 mm x 72 mm x 35 mm)
Weight:	4.6 oz (130 g) with batteries, 2.6 oz w/o batteries (73g)
Battery Type:	Two (2) AA non-rechargeable alkaline batteries (BAT 001), approx. 50 hours battery life; or Two (2) AA rechargeable NiMH batteries (BAT 026), 1500 mAh, approx. 32 hours battery life.
Operating Frequencies:	8 channels: 72.1, 72.3, 72.6, 72.8, 74.7, 75.4, 75.7, or 75.9 MHz.
Indicators:	Power on: bright green. Low Battery: flashes. Sleep Mode: fade on, fade off.
Receiver Antenna:	Integral with earphone/headphone cord
Output Connector:	3.5 mm stereo/mono headphone jack, accepts either stereo or mono plug
Approvals:	FCC, Industry Canada, RoHS, WEEE
Warranty:	Lifetime PLUS Limited Warranty, 90 Days on most accessories.

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Please do not dispose of equipment or used batteries in the household trash. Take the equipment to an electronics recycling center or return to the factory. Take batteries to a retail or community recycling location.

Battery Information

The PPA R37-8 uses 2 AA batteries. We recommend using Williams AV Alkaline batteries (BAT 001-2) or rechargeable NiMH batteries (BAT 026-2) for longer battery life than some other brand batteries.

Low Battery Indicator

The “On” indicator on top of the unit will flash to indicate low battery.

Multi Charger CHG 3512

This optional charger is ideal for customers using rechargeable batteries. It will charge up to 12 receivers at one time. See our website or call us for ordering information.

Accessories

To order accessories including earphones and replacement earpads, batteries, and chargers, contact Williams AV at **800-328-6190** or order on-line at: **www.williamsav.com**.

Lifetime PLUS Limited Warranty

Williams AV products are engineered, designed, and manufactured under carefully controlled conditions to provide you with many years of reliable service.

Williams AV warrants the PPA R37-8 against defects in materials and workmanship under normal use and conditions for the useful lifetime of the product from date of purchase.

This warranty is available to the original end purchaser of the product and CAN BE transferred to subsequent purchasers of the product.

Microphones, earphones, headphones, batteries, chargers, cables, carry cases, and most other accessory products carry a 90-day warranty.

Williams AV has no control over the conditions under which this product is used. Williams AV, therefore, disclaims all warranties not set forth above, both express and implied, with respect to the Personal PA[®] FM Listening System, including but not limited to, any implied warranty of merchantability or fitness of use of such equipment including, without limitation, any warranty that the use of such equipment for any purpose will comply with applicable laws and regulations. Williams AV shall not be liable to any person or entity for any medical expenses or any direct, incidental or consequential damages caused by any use, defect, failure or malfunctioning of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise, the sole remedy for any defect, failure or malfunction of the products is replacement of the product. No person has any authority to bind Williams AV to any representation or warranty with respect to the Personal PA[®] FM Listening System. Unauthorized repairs or modifications will void the warranty. This warranty is void if damage occurred because of misuse, or if the product has been repaired or modified by anyone other than a factory authorized service technician. Warranty does not cover normal wear and tear on the product or any other physical damage unless the damage was the result of a manufacturing defect. Williams AV is not liable for consequential damages due to any failure of equipment to perform as intended. Williams AV shall bear no responsibility or obligation with respect to the manner of use of any equipment sold by it.

This warranty does not cover reimbursement for your costs of removing and transporting the product for warranty service evaluation or installation of any replacement product provided under this warranty.

The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Limited Warranty did not contain the particular part or term held to be invalid. The terms of the warranty are governed by the laws of the State of Minnesota.

Prices and the specifications of the products are subject to change without notice.

For Complete Warranty Statement go to: www.williamsav.com/warranty-statement

Customer Service

If you experience difficulty with your system, call Toll-Free for Customer Assistance:

1-800-328-6190 (U.S.A.) or

+1 952 943 2252 (Outside the U.S.A.)

If it is necessary to return the system for service, your Customer Service Representative will give you a Return Authorization Number (RA) and shipping instructions.

Pack the system carefully and send it to:

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