Hear what you want to hear by amplifying the sounds and voices most important to you.

1

Plug your earphone or headphone into the headphone jack.

2

Position the microphone near the preferred sound.

3

Adjust the volume to your comfort.

Great for



One-on-one discussions



Television listening



Conversing in the car



Amplify your life with the Pocketalker®

Hear more, connect more with the Williams Sound Pocketalker personal amplifier.

Don't let your important moments go unheard. The Pocketalker ensures you never miss a single word, whether you're in a noisy environment or enjoying a quiet conversation.

Personally tailor your listening experiences with this out-of-the-box hearing assistance solution. Stay connected and involved by amplifying the sounds and voices most important to you - at a fraction of the cost of hearing aids.

WILLIAMSSOUND

A Williams AV Brand

Hear the **Important** things in life



1.800.328.6190 (USA) | +1.952.943.2252 (World)

info@williamsav.com

www.williamsav.com/personal

Copyright 2023 - Williams AV LLC BROM073A

P O C K E T A L K E R[®]

Personal Amplifier



66

I bought the Pocketalker for my
89 year old father and was in
tears at the look on his face when
he could hear and understand
me again. We put on some old
music he used to love and he just
sat there with a smile on his face.
Thank you for giving my father
back to me!

Debbie M

Verified Customer Review



Pocketalker vs. Hearing Aids

A hearing aid is used to amplify all sounds. It is prescribed and tailored based on an individual's audiogram and may require periodic adjustments by an audiologist.

The **Pocketalker** personal amplifier is used to **amplify specific sounds or voices**.

Simple controls allow users to adjust tone and balance to their needs without requiring a prescription or help from an audiologist.



Hear the important things in life

Simple Design. **Personalized Listening.**

The Pocketalker is designed to provide a personalized, easy-to-use listening experience.

Use the control to adjust the volume of the sound up to 120 decibels, which is twice the level of a normal conversation.

Further personalize your listening experience to enhance sound clarity using tone control to adjust low and high-frequency sounds and the balance control to adjust the sound levels in each ear.



Easy-to-use intuitive design



Long-lasting power



Adjustable tone control for improved clarity



Ergonomic, lightweight design



Simple maintenance

Pocketalker Comparison Chart



	Ultra	2.0
5-Year Warranty	J	J
Tone Control	J	J
Balance Adjustment		J
Built-in Telecoil (for hearing looped venues)		J
Rechargeable Batteries (Optional)		Up to 80 hr. Runtime
AAA Alkaline Batteries	200 hr. Runtime	105 hr. Runtime
TV Extension Cable	J	J
Internal Microphone (works without plug-in microphone)		J
Accessories	Single-ear headphones and earbud Lanyard neck	Dual-ear headphones and earbuds Lanyard neck
	strap with safety clasp	strap with safety clasp